Thoraflex™ Hybrid Device

Instructions for Use

CE
0086

 **VASCUTEK**
TERUMO

English	Instructions for Use	4
Français	Mode d'emploi	24
Deutsch	Gebrauchsanweisung	44
Nederlands	Gebruiksaanwijzing	64
Italiano	Istruzioni per l'uso	84
Español	Instrucciones de uso	104
Português	Instruções de Utilização	124
Dansk	Brugervejledning	144
Norsk	Bruksanvisning	164
Česky	Návod k použití	184
Slovensky	Návod na použitie	204
Русский	Инструкция по эксплуатации	224

English

Vascutek Ltd. Thoraflex™ Hybrid Device

Instructions for Use

CONTENTS

SECTION 1 INSTRUCTIONS FOR USE

- 1.1 Components and Description
- 1.2 Indications
- 1.3 Contraindications
- 1.4 Cautions
- 1.5 Magnetic Resonance Imaging (MRI) Safety

SECTION 2 SIZING

SECTION 3 INSTRUCTIONS FOR IMPLANTATION OF THE THORAFLEX™ HYBRID DEVICE

- 3.1 Preparation of the Thoraflex™ Hybrid Delivery System
- 3.2 Forming of the Thoraflex™ Hybrid Delivery System
- 3.3 Introduction of the Thoraflex™ Hybrid Delivery System Using a Guidewire
- 3.4 Positioning the Thoraflex™ Hybrid Delivery System
- 3.5 Deployment Sequence for the Thoraflex™ Hybrid Device
- 3.6 Thoraflex™ Hybrid Device Anastomoses

SECTION 4 ADDITIONAL INFORMATION

- 4.1 Clinical Experience Summary
- 4.2 Use of the Thoraflex™ Hybrid Device Perfusion Side Branch
- 4.3 Origin of Gelatin
- 4.4 Sterilisation
- 4.5 Packaging
- 4.6 Additional Labels
- 4.7 Disposal of the Thoraflex™ Hybrid Delivery System
- 4.8 Returning a Thoraflex™ Hybrid Delivery System or Device

SECTION 5 TROUBLESHOOTING GUIDE

- 5.1 Difficulty Advancing the Delivery System to the Intended Deployment Site
- 5.2 Sheath Retraction - A High Force is Required to Retract the Sheath and Deploy the Device
- 5.3 Sheath Retraction – The Strap becomes Detached from the Delivery Sheath
- 5.4 Sheath Retraction – The Strap Handle Breaks in Two
- 5.5 The Splitter is Detached During / Prior to Deployment
- 5.6 The Handle Release Clip is not Attached to the Delivery System Handle
- 5.7 Difficulty Removing the Guidewire
- 5.8 Difficulty Removing the Release Wire from the Delivery System
- 5.9 Difficulty Removing the Delivery System – High Angulation
- 5.10 Difficulty Removing the Delivery System– Release Wire and Guidewire

SECTION 1 INSTRUCTIONS FOR USE

This booklet provides instructions for the use of the Thoraflex™ Hybrid System.

SECTION 1.1 COMPONENTS AND DESCRIPTION

The Vascutek Ltd. Thoraflex™ Hybrid System consists of the Thoraflex™ Hybrid Device contained within a delivery system as shown in Figure 1. The Thoraflex™ Hybrid Device is made from gelatin sealed woven polyester, a section of which has the addition of nitinol ring stents as shown in Figures 2 and 3. The device is specifically designed to repair or replace damaged or diseased vessels of the aortic arch and descending aorta with or without involvement of the ascending aorta in cases of aneurysm and/or dissection by open surgical repair. The branches of the Plexus 4 version (Figure 2) are provided to accommodate reconstruction of the major aortic branch vessels and the intra-operative attachment of a perfusion cannula during cardio-pulmonary bypass where antegrade perfusion techniques are used. The stented section of the device enables treatment of the diseased descending thoracic aorta without the need to anastomose a stented and unstented device together.

The Thoraflex™ Hybrid Device has been impregnated with an absorbable protein. The aim of the impregnation is to provide a polyester vascular device that does not require pre-clotting. The protein is a bovine gelatin which has been cross-linked to a set level in order to control its rate of removal. It serves in place of fibrin, which seals the polyester device during normal pre-clotting. The gelatin is hydrolysed within 14 days and is replaced by normal tissue incorporation. Gelatin was chosen as it is a non-toxic protein, a fact which is reflected by its extensive use as a safe plasma expander.


Figure 1 – Thoraflex™ Hybrid Delivery System (Plexus 4 Version)


Figure 2 – Thoraflex™ Hybrid Plexus 4 Device


Figure 3 – Thoraflex™ Hybrid Ante-Flo Device

SECTION 1.2 INDICATIONS

The Thoraflex™ Hybrid Device is intended for single stage open surgical repair or replacement of damaged or diseased vessels of the aortic arch and descending aorta with or without involvement of the ascending aorta in cases of aneurysm and/or dissection.

SECTION 1.3 CONTRAINDICATIONS

This device should not be implanted in patients who exhibit;

1. Sensitivity to Polyester, Nitinol, or materials of bovine origin.
2. Active infection.

SECTION 1.4 CAUTIONS

1. The vascular graft material is based on a woven structure and therefore must be cut with a cautery to minimise fraying.

NOTE: IMMERSION OF THE DEVICE IN A SALINE SOLUTION IMMEDIATELY PRIOR TO USE WILL PREVENT FOCAL BURNING WHICH MAY RESULT DURING CAUTERISATION. The device must be immersed in a saline solution for at least 1 minute, but no longer than 5 minutes. The device must not be allowed to dry out after soaking.

2. It is of particular importance that the stented section is pre-soaked as advised in section 3.1 as this will greatly reduce the force required to deploy the device.

3. **DO NOT PRECLOT.** These devices contain a gelatin sealant and must not be pre-clotted.

4. DO NOT USE BEYOND THE INDICATED EXPIRATION DATE. The gelatin impregnation may not meet the design specification after the expiration date.

5. DO NOT RESTERILISE. FOR SINGLE USE ONLY. Do not reuse, reprocess or re-sterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in deterioration of health or death of patients. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross infection, including but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient end-user.

6. Store in clean, dry area at room temperature.

7. The device must be implanted within one month after removal from the foil pouch.

8. Clamping may damage the prosthesis. Atraumatic clamps, ideally with soft shod jaws, should be used with a minimum application of force. Excessive force should be avoided as it will damage the polyester fibres and the gelatin impregnation.

9. Excessive tension on the prosthesis should be avoided.

10. Round body taper point needles should be used when implanting the device to minimise fibre damage.

11. If de-airing is required then the smallest needle possible should be used. A 19 gauge needle is normally sufficient. Hypodermic needles have a cutting point which may result in blood leakage and require repair by suturing.

12. Devices should be selected according to the Vascutek Ltd. Thoraflex™ Hybrid Device Sizing Chart (see section 2).

13. Long term performance of the graft has not been established; therefore patients should be monitored on a regular basis for adverse events, for example endoleaks and aneurysm growth.

14. The use of a balloon expandable stent e.g. Palmaz® stent, to treat an endoleak may result in abrasion of the graft material leading to graft failure or fatigue.

15. Excess delivery system angulation will cause more kinking of the sheath and therefore require an increased deployment force.

16. Delivery system removal - if the system is introduced around a curve it must be removed following an identical path in order to avoid moving the device or causing trauma to the vessel.

17. In dissection cases, additional care must be taken during the insertion and removal of the delivery system in order to minimise the risk of trauma to the vessel wall.

18. The graft section of the Thoraflex™ Hybrid Device is made from gelatin sealed Gelweave™ polyester material. Due to its use in conjunction with a delivery system, the graft section may suffer some slight initial blood loss when compared to a standard gelatin sealed Gelweave™ graft.

19. The manufacturing process for gelatin sealed vascular grafts uses the cross-linking agent formaldehyde to achieve the graft performance. All gelatin sealed grafts are thoroughly rinsed with RO water to reduce residual formaldehyde, however, residual amounts may be present in the finished graft. Formaldehyde is also found at low levels naturally in the body, some of which is derived from food. Formaldehyde is known to be mutagenic and carcinogenic. The risks of these potential harms from the product have not been established clinically.

20. The Thoraflex™ Hybrid Device has not been tested in paediatric populations.

21. IF A GUIDEWIRE WAS USED DURING THE DEPLOYMENT OF THE DEVICE IT MUST BE REMOVED FROM THE DELIVERY SYSTEM BEFORE THE RELEASE WIRE IS REMOVED – See Section 3.5.3.

SECTION 1.5 MAGNETIC RESONANCE IMAGING (MRI) SAFETY

The Thoraflex™ Hybrid Device was determined to be **Magnetic Resonance (MR)-conditional**. Non-clinical testing determined that grafts with radiopaque markers were MR conditional. A patient with this device can be scanned safely, immediately after placement of the device under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla or less.
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less.

MAGNETIC RESONANCE IMAGING (MRI) RELATED HEATING

In non-clinical testing, grafts with radiopaque markers produced the following temperature rises during Magnetic Resonance Imaging performed for 15 minutes of scanning (per pulse sequence) in 1.5 Tesla/64 MHz (Magnetom, Siemens Medical Solutions, Malvern PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	1,5 Tesla	3 Tesla
MR system reported, whole body averaged SAR	2,9 W/kg	2,9 W/kg
Calorimetry measured values, whole body averaged SAR	2,1 W/kg	2,7 W/kg
Highest temperature change	+1,7°C	+2,0°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artefact Information

MR image quality may be compromised if the area of interest is in the exact same area or close to the position of the graft with radiopaque markers. Therefore, optimisation of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artefact size (i.e. as seen on the gradient echo pulse sequence) extends approximately 10mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	15, 818 mm ²	1, 424 mm ²	19, 077 mm ²	2, 012 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

This information is based on information from the Food and Drug Administration and the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marketing Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

SECTION 2 SIZING

The Thoraflex™ Hybrid Sizing Chart incorporates a 15-25% oversize of ring stent diameter to aortic diameter. Aortic diameter is based on inner vessel diameter (ID) measurements therefore no further oversize is required. If outside vessel diameters (OD) are measured then an allowance for the vessel wall diameter must be made before using the sizing chart for device selection.

Vascutek has performed testing and recommends that a 15-25% oversize and a 40mm distal landing zone will provide optimum sealing within healthy vessel (Figure 4a, 4b, 4c and 5).

When oversizing in a dissection clinical judgement must be used on an individual patient basis.

NOTE: The diameter of the sheath measures 10mm.

Catalogue No. (Plexus Design)	Catalogue No. (Ante-Flo Design)	D1 Graft ID (mm)	D2 Stent Graft OD (mm)	D3* Descending Landing Zone Vessel ID (mm)	Branch Configuration IA, LCC, LSA (mm)	L1** Stent Graft Nominal Length (mm)
THP2224X100B	THA2224X100B	22	24	19 - 21	10, 8, 8	100
THP2426X100B	THA2426X100B	24	26	20 - 22	10, 8, 8	100
THP2628X100B	THA2628X100B	26	28	22 - 24	12, 8, 10	100
THP2830X100B	THA2830X100B	28	30	24 - 26	12, 8, 10	100
THP3032X100B	THA3032X100B	30	32	25 - 27	12, 8, 10	100
THP3034X100B	THA3034X100B	30	34	27 - 29	12, 8, 10	100
THP3036X100B	THA3036X100B	30	36	29 - 31	12, 8, 10	100
THP3038X100B	THA3038X100B	30	38	30 - 33	12, 8, 10	100
THP3040X100B	THA3040X100B	30	40	32 - 34	12, 8, 10	100
THP3240X100B	THA3240X100B	32	40	32 - 34	12, 8, 10	100
THP2224X150B	THA2224X150B	22	24	19 - 21	10, 8, 8	150
THP2426X150B	THA2426X150B	24	26	20 - 22	10, 8, 8	150
THP2628X150B	THA2628X150B	26	28	22 - 24	12, 8, 10	150
THP2830X150B	THA2830X150B	28	30	24 - 26	12, 8, 10	150
THP3032X150B	THA3032X150B	30	32	25 - 27	12, 8, 10	150
THP3034X150B	THA3034X150B	30	34	27 - 29	12, 8, 10	150
THP3036X150B	THA3036X150B	30	36	29 - 31	12, 8, 10	150
THP3038X150B	THA3038X150B	30	38	30 - 33	12, 8, 10	150
THP3040X150B	THA3040X150B	30	40	32 - 34	12, 8, 10	150
THP3240X150B	THA3240X150B	32	40	32 - 34	12, 8, 10	150

Figure 4a - Thoraflex™ Hybrid Device Sizing Chart

*See Figure 5. Oversize stated is approximately 15-25%.

**Nominal length quoted.

It should be considered when measuring length that the device will take the outer curve of the aorta. Refer to Figure 5.

Some movement of the distal ring of the Thoraflex™ Hybrid Device may occur following re-perfusion of the thoracic aorta.


Figure 4b – Thoraflex™ Hybrid Plexus 4 Device Dimensions

***Length shown to indicate marker position. Ring configuration will vary.


Figure 4c – Thoraflex™ Hybrid Ante-Flo Device Dimensions

***Length shown to indicate marker position. Ring configuration will vary.


Figure 5 – Device Sizing

Some movement of the distal ring of the Thoraflex™ Hybrid Device may occur following re-perfusion of the thoracic aorta.

SECTION 3 INSTRUCTIONS FOR IMPLANTATION OF THE THORAFLEX™ HYBRID DEVICE

SECTION 3.1 PREPARATION OF THE THORAFLEX™ HYBRID DELIVERY SYSTEM

Once in the sterile field, the entire Thoraflex™ Hybrid Delivery System must be thoroughly pre-soaked for at least 1 minute but no longer than 5 minutes in a saline solution to ensure adequate absorption.

Note: The delivery system and the device must not be allowed to dry out after soaking. Pre-soaking will reduce the force necessary to unsheathe the stented section of the device.

SECTION 3.2 FORMING OF THE THORAFLEX™ HYBRID DELIVERY SYSTEM

The stented section of the Thoraflex™ Hybrid Delivery System (Figure 6) can be shaped to resemble the anatomy of the aorta in the area of the stent graft only.

CAUTION: DO NOT BEND THE STENTED SECTION WITHIN 10MM OF THE SPLITTER OR WHILE HOLDING THE HANDLE OF THE SYSTEM.

CAUTION: IF THE SYSTEM IS FORMED INTO AN ANGLE OF GREATER THAN 50° A HIGHER FORCE WILL BE NEEDED TO DEPLOY THE STENT.


Figure 6 – Shaping of the Stented Section of the System

If any major sheath kinking is evident after shaping the delivery system then local pressure should be applied to the kinks in order to reduce the sheath folding and remove any sharp points (Figure 7).

Figure 7 – Removal of Kinks in the Sheath

SECTION 3.3 INTRODUCTION OF THE THORAFLEX™ HYBRID DELIVERY SYSTEM USING A GUIDEWIRES

It is recommended that the Thoraflex™ Hybrid Delivery System is used with a guidewire (Figure 8). The tip contains a choice of two guidewire access ports (Figure 9). The guidewire can be fed through either port, and then along the outside of the sheath. The delivery system can then be moved along it into position.


Figure 8 – Guidewire Use in the Aorta


Figure 9 – Insertion of the Guidewire Through Tip Guidewire Access Ports

SECTION 3.4 POSITIONING THE THORAFLEX™ HYBRID DELIVERY SYSTEM

The Thoraflex™ Hybrid Delivery System must be placed through the opened aortic arch into the descending thoracic aorta. This should be done over a guidewire in order to ensure that the correct lumen is being treated, for example in cases of dissection (Figure 10). When positioning the delivery system ensure that the splitter release clip is accessible and the collar is positioned correctly in relation to the anastomotic site.

The splitter should be positioned in the distal aorta so that when the device is deployed the collar is in the correct position (Figure 10). For the Thoraflex™ Hybrid Plexus 4 version, the Delivery System should be orientated so that the device branches and aortic arch vessels are aligned.


Figure 10 – Positioning of the Thoraflex™ Hybrid Delivery System

SECTION 3.5 DEPLOYMENT SEQUENCE FOR THE THORAFLEX™ HYBRID DEVICE

3.5.1 Sheath Retraction (Device Release Stage I)

When the optimum orientation and position has been achieved the delivery system should be unsheathed. In order to unsheathe the device, firmly stabilise the handle with one hand and with the other hand pull back the strap in-line with the handle (Figure 11a, b, c). This will simultaneously retract and split the sheath allowing it to be completely removed from the delivery system. The complete stented section of the device will now be unsheathed.


Figure 11a –Sheath Retraction


Figure 11b –Sheath Retraction


Figure 11c –Sheath Retraction

3.5.2 Removal of the Sheath Splitter (Device Release Stage II)

Once the sheath has been removed, the splitter is detached from the delivery system by cutting the suture (Figure 12a). Ensure the graft fabric under the splitter is opened up to facilitate the removal of the handle (Figure 12b).


Figure 12a – Remove Splitter by Cutting the Suture

Figure 12b – Graft Fabric Opened Around the Collar

3.5.3 Guidewire Removal

CAUTION: IF A GUIDEWIRE WAS USED DURING THE DEPLOYMENT OF THE DEVICE IT MUST BE REMOVED FROM THE SYSTEM BEFORE THE RELEASE WIRE IS REMOVED (Figure 13). This step enables the guidewire to be removed while the device remains held in position by the system thereby preventing movement of the stented section.


Figure 13 – Guidewire Removal

3.5.4 Release Wire Removal (Device Release Stage III)

In order to fully release the device from the delivery system pull the red release clip and attached wire out of the delivery system handle (Figure 14). The release wire should be pulled out proximally, in line with the delivery system handle. The distal end of the stent graft will now be released from the delivery system (Figure 14).


Figure 14 – Release Wire Removal

3.5.5 Delivery System Removal

Once the device has been released from the delivery system the remaining handle assembly must be removed from the device. This is removed by gently pulling the handle proximally ensuring that the device is sufficiently loose around the shaft to allow removal without disturbing the graft. If the delivery system was introduced around a curve, as shown in Figure 15, it must be removed following the identical path in order to avoid moving the device or causing trauma to the vessel.


Figure 15 – Delivery System Removal

SECTION 3.6 THORAFLEX™ HYBRID DEVICE ANASTOMOSES

Once the delivery system has been removed, the collar should be sutured to the native aortic vessel in order to provide fixation and stability to the device. The exact technique used is at the discretion of the surgeon implanting the device however a circumferential anastomosis is required to ensure that the implant is sealed correctly (Figure 16). The remaining anastomoses are now carried out.

Note: Some movement of the distal ring of the Thoraflex™ Hybrid Device may occur following re-perfusion of the thoracic aorta.


Figure 16 – Thoraflex™ Hybrid Plexus 4 Device Anastomoses

SECTION 4 ADDITIONAL INFORMATION

SECTION 4.1 CLINICAL EXPERIENCE SUMMARY

Shrestha et al. (2016) published results of the first 100 patients treated with the Thoraflex™ Hybrid device at a single centre in Hannover, Germany. This cohort includes data from 34 patients treated as part of the initial pivotal pre-market study of the device undertaken at the same institution. There were 37 Acute dissection patients, 31 chronic dissection and 32 thoracic aortic aneurysm cases. Twenty-eight patients had undergone previous surgery (28%), whilst twelve (12%) presented with Marfan syndrome.

Postoperatively, re-thoracotomy for bleeding was required in ten patients (10%), whilst stroke was reported in nine (9%). The stroke rate was higher in acute aortic dissection patients (14%, n = 5) compared with chronic aortic dissection (10%, n = 3) and aneurysm patients (3%, n = 1). A total of 14 patients (14%) experienced acute renal failure requiring dialysis, and eight patients (8%) were discharged on dialysis. Paraparesis occurred in seven patients (7%), three were acute dissection patients who were not routinely treated with cerebrospinal fluid (CSF) drainage. All but one recovered with CSF drainage post-operatively. The perioperative mortality rate, including 30-day mortality and in-hospital deaths, was 7% (seven patients).

Follow up was 100% complete. During a mean follow up of 2.7 years (range 1.8 to 4.6 years), 22 patients required further procedures in the downstream aorta. Of these 22 procedures, 15 were anticipated/planned extensions (68%). 13 of the extension procedures were performed using an endovascular approach (59%). Death following discharge was reported in 12 patients (13%) including three acute dissection patients (9%), four chronic dissection patients (13%) and five aneurysm patients (17%).

SHRESTHA, M., KAUFELD, T., BECKMANN, E., FLEISSNER, F., UMMINGER, J., ABD ALHADI, F., BOETHIG, D., KRUEGER, H., HAVERICH, A. & MARTENS, A. 2016. Total aortic arch replacement with a novel 4-branched frozen elephant trunk prosthesis: Single-center results of the first 100 patients. *J Thorac Cardiovasc Surg*, 152, 148-159 e1.

Based on data from the published literature and the pivotal trial of the Thoraflex Hybrid device, Vascutek have identified the following potential procedural risks/complications which can be associated with Frozen Elephant Trunk surgery:

Perioperative/In-Hospital Mortality, Rethoracotomy for Bleeding, Stroke, TIA, SCI/Paraparesis/Paraplegia, Recurrent Nerve Palsy, Prolonged Ventilation, Respiratory Complications, Cardiac Complications, Acute Kidney Injury, Renal Failure/Dialysis, Requirement for Extension of Thoraflex Hybrid Device (via open or endovascular procedure), Endoleak, Aortic Rupture, Pulmonary Embolism, Systemic/Deep Infection, Wound Infection

SECTION 4.2 USE OF THE THORAFLEX™ HYBRID DEVICE PERfusion SIDE BRANCH

Initiation of Antegrade Perfusion: The bypass catheter should be placed in the side branch of the Plexus 4 or Ante-Flo device and securely attached.

Completion of Antegrade Perfusion: Once bypass is complete, the cannula side branch of the Plexus 4 or Ante-Flo device should be cut off and the remaining stump over-sewn using standard surgical technique.

SECTION 4.3 ORIGIN OF GELATIN

Vascutek uses gelatin manufactured from animals native to, and exclusively raised in the United States of America. The United States of America is classified as a negligible BSE risk country according to the OIE categorisation (adopted by the European Union with the Regulation (EC) N°722/2007). The gelatin is hydrolysed within approximately 14 days and is replaced by normal tissue incorporation.

SECTION 4.4 STERILISATION

These systems have been sterilised using Ethylene Oxide and are supplied sterile. The Tyvek® seal on both intermediate and inner pouches must be intact. Any damage to the pouches renders the system non-sterile.

Note: In the event of damage to the primary packaging, the product must not be used and should be returned immediately to the supplier.

CAUTION: THE THORAFLEX™ HYBRID SYSTEM MUST NOT BE RE-STERILISED.

SECTION 4.5 PACKAGING

Tyvek® pouches are enclosed in a foil pouch that serves as a vapour barrier and preserves optimal system characteristics. A sachet containing a desiccant is included to aid this purpose.

Note: The foil pouch and outer Tyvek® pouch are not sterile. Only the innermost pouch and tray can be introduced to the sterile field.

SECTION 4.6 ADDITIONAL LABELS

Additional labels are enclosed for use on patient records.

SECTION 4.7 DISPOSAL OF THE THORAFLEX™ HYBRID DELIVERY SYSTEM

At the end of the procedure care must be taken to ensure safe disposal of the Thoraflex™ Hybrid Delivery System. Each operating team must ensure local and national regulatory requirements for the disposal of contaminated clinical waste products are adhered to.

SECTION 4.8 RETURNING A THORAFLEX™ HYBRID DELIVERY SYSTEM OR DEVICE


All explanted delivery systems/devices should be returned to Vascutek for analysis as soon as possible. In the event of a used delivery system/device needing to be returned to Vascutek, it is a requirement to have the delivery system/device, and any other items used in the procedure to be returned in an explants box which can be obtained from Vascutek's Quality Assurance Department. If required, explant kits can be requested at complaints@vascutek.com or through your local distributor and will be provided for the retrieval and preservation of the explanted device and/or delivery system or other components for transit to Vascutek.

SECTION 5 TROUBLESHOOTING GUIDE

5.1 DIFFICULTY ADVANCING THE DELIVERY SYSTEM TO THE INTENDED DEPLOYMENT SITE

Potential Causes

- A. The level of angulation of the sheath/malleable shaft may not be compatible with the aortic anatomy.
- B. During introduction the system does not follow the path of the guidewire.


Recommended Actions

- A.1. Remove the delivery system, check the degree of angulation of the delivery system shaft is suitable for the anatomy.

- B.1. Feed the guidewire through the other tip guidewire port. This should minimise the effect of this issue.

Note: The guidewire does not travel through the centre of the delivery system. As a result consideration must be taken of the appropriate angle when advancing the system.

- B.2. If A1 and B1 do not work then use another device.

5.2 SHEATH RETRACTION - A HIGH FORCE IS REQUIRED TO RETRACT THE SHEATH AND DEPLOY THE DEVICE


Potential Cause

High level of angulation of the sheath and malleable shaft may result in excessive kinking of the sheath. This can lead to higher deployment forces being experienced during sheath retraction.

Recommended Actions

Where high deployment forces are experienced, applying a syringe like grip to the system can provide greater control for the initial stage of the retraction process.

The following steps may also resolve this issue.

1. Remove the system and inspect the sheath for any excessive kinking.
2. Where possible straighten the delivery system to remove the kinks.
3. Kinks can also be removed from the sheath by using localised compression of the affected area.
4. Repeat steps 1-3 until the deployment force required has been reduced.


Note: Ensure one hand stabilises the handle. The graft section should be moved forward thus providing a longer length of handle to hold.

5.3 SHEATH RETRACTION - THE STRAP BECOMES DETACHED FROM THE DELIVERY SHEATH

Potential Cause


Excessively high deployment forces or poor strap to sheath attachment strength may result in the detachment of the strap from the sheath.

Recommended Actions

1. Attach one pair of atraumatic artery forceps to each side of the sheath.
2. Ensure the handle is stabilised.
3. With the handle stabilised, pull the forceps to retract both parts of the sheath simultaneously.


Note: Both parts of the sheath must be pulled back simultaneously to ensure the correct unsheathing of the stented section of the device.

5.4 SHEATH RETRACTION - THE STRAP HANDLE BREAKS IN TWO

Potential Cause


An excessive deployment force can contribute to the potential fracture and/or weakening of the strap. This may cause it to break during the device deployment.

Recommended Actions

1. Gather both halves of the strap in one hand.
2. Stabilise the handle of the delivery system.
3. Retract both halves of strap together thus splitting the sheath in two.


Note: This process could be completed using atraumatic artery forceps as described in section 5.3 Sheath Retraction - The strap becomes detached from the delivery sheath

5.5 THE SPLITTER IS DETACHED DURING / PRIOR TO DEPLOYMENT

Potential Cause

Splitter suture is loose or unfastened during use.


Recommended Actions

1. The splitter can be manually closed and held during deployment or if partially open can be reattached.
2. Ensure that the fabric of the graft is not caught in the splitter when re-attaching.
3. Push both sides of the splitter together until they click together.
4. The splitter should be held closed during the sheath retraction element of the procedure.

Note: If the hinge of the splitter breaks care must be taken to ensure that both parts of the splitter are removed from the chest cavity after use.

Note: The splitter is not essential to deploy the device.

5.6 THE HANDLE RELEASE CLIP IS NOT ATTACHED TO THE DELIVERY SYSTEM HANDLE

Potential Causes

- A. The red release clip has become loose, but is still attached to the handle, either during transit or use.
- B. The red release clip is no longer attached to the release wire.

Recommended Actions

- A.1. Check that the red release clip is still attached to its release wire by gently pulling it to ensure that tension can still be felt.
- A.2. If still attached, the red release clip can be pushed back into its housing in the handle.
- A.3. The red release clip can be removed and retracted as per the IFU section 3.5.4.

- B.1. Remove the release wire

- B.2. The release wire should be retracted as described in the IFU section 3.5.4 usingatraumatic artery forceps. Theatraumatic artery forceps should be attached to the release wire and retracted.

Note: Only one wire is present in the delivery system

5.7 DIFFICULTY REMOVING THE GUIDEWIRE

Potential Cause


The combination of delivery system manipulation and the guidewire port used has prevented the easy removal of the guidewire from the delivery system.

Recommended Actions

1. Ensure that the handle is stabilised during guidewire removal to avoid compromising the position and stability of the stented section of the device.
2. Try to remove the guidewire in the opposite direction to how it was inserted i.e. distal removal not proximal / proximal removal not distal.


Note: The direction of guidewire removal may be compromised by the anatomy.

5.8 DIFFICULTY REMOVING THE RELEASE WIRE FROM THE DELIVERY SYSTEM

Potential Cause

If resistance is felt by the operator when removing the release wire, there is the possibility that the release wire has become kinked and/or snapped within the device. The snapped portion of the release wire must be removed before the deployment can be completed.

Recommended Actions

1. The release wire should be inspected once it has been removed from the system. Look for any sign of fracture, damage or shortened length.


Final length should be approximately;
• Cat No. TH.....X100B = 530mm
• Cat No. TH.....X150B = 630mm

2. Check the wire entry ports on the handle of the delivery system. If a piece of wire is still there remove it using forceps.

The release wire runs from the release clip to the tip of the delivery system and terminates back inside the handle.


If no wire is present on either side of the handle then the device should have been freed from the delivery system. Any difficulty experienced must be due to the angulation of the system. See section 5.9

5.9 DIFFICULTY REMOVING THE DELIVERY SYSTEM – HIGH ANGULATION

Potential Cause

The level of angulation applied to the delivery system is causing difficulty during handle removal.

Recommended Actions

1. Take account of the curvature applied to the delivery system before application, the removal may have to follow a similar path.
2. Remove gently.


5.10 DIFFICULTY REMOVING THE DELIVERY SYSTEM – RELEASE WIRE AND GUIDEWIRE


Potential Causes

- A. The delivery system is still attached to the device.
- B. The delivery system is still attached to the guidewire.

Note: If the handle is removed whilst still attached to the guidewire, the stented section of the device will become dislodged


- A.1. Check that the release wire and splitter have been removed.
- A.2. Ensure the graft has been unfolded in the area where the splitter was.
- A.3. Gently move the handle distally before removing.

- B.1. If the device is positioned correctly (sheathed or unsheathed) and is still connected to the delivery system, the guidewire can be removed without dislodging the device

Note: This is the last point in the procedure where the guidewire can be removed without significant consequence or detrimental effect to the device.


5.10 DIFFICULTY REMOVING THE DELIVERY SYSTEM – RELEASE WIRE AND GUIDEWIRE

B.2. After removal of the release wire, the device is no longer attached to the delivery system, hence removal of the guidewire may dislodge the device if not previously removed.

B.3. If the handle is removed while attached to the guidewire, this may cause the device to be moved from the deployment site. The effects may be detrimental without secondary intervention.

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VASCUTEK Ltd