REF

T401231, T401231LS, T401231U, T401232, T401232U

Titan™ Stabilizer Attachments

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Titan™ Stabilizör Eklentileri

This page is intentionally left blank.

Titan™ Stabilizer Attachments

Table of Contents

English	5
Česky	11
Dansk	17
Nederlands	23
Suomi	29
Français (France)	35
Deutsch	41
Magyar	47
Italiano	53
Norsk	59
Polski	65
Português (Brasil)	71
Русский язык	77
简体中文	83
Slovenčina	89
Español (España)	95
Svenska	101
繁體中文	107
Türkçe	113
Graphical Symbols & Definitions	118

This page is intentionally left blank.

English

Instructions for Use

Titan™ Stabilizer Attachments

Catalogue Numbers T401231, T401231LS, T401231U, T401232, T401232U

Compatible with Terumo Stabilizer Arms: T401160, T401161, T401162, T401163

CAUTION: Before use, please read carefully the instructions for use of the Terumo Stabilizer Arms in addition to these instructions for use.

PRODUCT DESCRIPTION

The Titan™ Stabilizers are sterile single-use accessory devices for use with the Terumo Stabilizing Arms. The device provides stabilization and coronary artery isolation. The Titan Stabilizer Attachments include a vacuum option to maintain the position of the device.

Included in the sterile package is:

Titan consisting of the following parts (see Figure 1):

- The stabilizer attachment assembly with an approximately 24" length of $\frac{1}{4}$ " outer diameter tubing with female luer connector on one end to connect to stopcock.
- Hi-flow 4-way stopcock

Accessory components consisting of the following parts:

- A 12' length of $\frac{1}{4}$ " internal diameter tubing with a male luer connector on one end and vacuum source connector on the other end.
- A 3' segment of $\frac{1}{4}$ " internal diameter tubing with female connectors on both sides to connect the canister to the vacuum regulator if required.
- Tubing Clip


Figure 1 – Titan™ Stabilizer Attachment with tubing and stopcock

INTENDED USE

The stabilizers are intended for use in beating heart CABG (coronary artery bypass graft) procedures to stabilize the anastomotic site (coronary artery isolation).

CONTRAINDICATIONS

The stabilizers are contraindicated for use in patients with local or systemic infection.

Do not position the stabilizers on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterile and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- Do not use product if packaging or product is damaged as sterility or integrity of the device may be compromised.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- **Use of excessive negative vacuum pressure may cause tissue damage.**

CAUTION

United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

1. Inspect the packaging and product thoroughly for damage prior to use. Using aseptic technique, remove the contents and pass to the surgical field for attachment to the compatible Terumo Stabilizer Arm.
- For connection to the Terumo Stabilizer Arm shaft configuration: Slide the sleeve of the arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with the logo on the sleeve of the Terumo Stabilizer Arm (see Figure 2). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 2 – Attachment of the Titan™ Stabilizer Attachment to the Stabilizer Arm

- For connection to the Terumo Stabilizer Arm ball configuration: Grasp the sleeve of the arm and press the ball of the attachment into the collet of the Stabilizer Arm (see Figure 3). Check for secure connection.

NOTE: Ensure that the sleeve is grasped during attachment insertion. Failure to grasp the sleeve may result in an inability to insert the attachment.


Figure 3 – Attachment of the Titan™ 360 Stabilizer Attachment to the Stabilizer Arm ball configuration

2. The shaft and U-Tube malleable zone can be bent to conform to anatomic requirements.
 - The shaft malleable zone (see Figure 4) can be bent to a maximum of 45° upwards and 45° downwards. Exceeding this may cause the shaft to break. Visually verify the integrity of the device after bending.
 - The U-Tube malleable zone (see Figure 4) can be bent upwards to a maximum of 45° and downwards to a maximum of 20°. Exceeding this may cause the U-Tube to break. Visually verify the integrity of the device after bending.
3. For the Titan™ Flex Stabilizer (T401232, T401232U), the U-Tube malleable zone can be bent to conform to anatomic requirements.
 - The Titan™ Flex Stabilizer's U-Tube malleable zone (see Figure 4) can be opened symmetrically to a maximum of 60° in the outwards direction. Exceeding this may cause the U-Tube to break. Visually verify the integrity of the device after bending.
4. For the Titan™ Flex Stabilizer (T401232, T401232U), the Foot malleable zone can be bent to conform to anatomic requirements.
 - The Foot malleable zone (Titan™ Flex only, see Figure 4) can be bent to a maximum of 20° downwards. Exceeding this may cause the Foot to break or device occlusion. Visually verify the integrity of the device after bending.


Figure 4 – Malleable Zones on the Titan™ Flex Stabilizer Attachments

WARNING

- Repeated bending and/or excessive flexing of the malleable zones will eventually cause the device to break. Slowly bend the malleable zones during manipulation. The integrity of the malleable zones should be visually verified after bending.
 - Take care not to occlude the vacuum lumen when bending devices.
5. Position the stabilizer attachment as required. Turn the main handle of the Terumo Stabilizing Arm clockwise to fasten the arm. Refer to the instructions for use of the Terumo Stabilizer Arms.

Use of vacuum option with Titan Stabilizer Attachments:

1. Attach male luer connection of accessory vacuum tubing to the female luer port of the vacuum control stopcock, making sure the vacuum source is switched OFF. Attach the vacuum source connector of the accessory vacuum tubing to the vacuum source; turn the vacuum control stopcock OFF to the vacuum source.
2. Turn the vacuum regulator of the vacuum source on to sufficient negative pressure to affix the Titan Stabilizer Attachment to tissue. Due to variations in OR equipment and utilities, use the least amount of vacuum necessary for desired product performance. The regulator should not exceed -400 mmHg.
3. Turn the vacuum control stopcock OPEN to the vacuum source and attachment once proper position has been obtained. Verify attachment of the stabilizer to the tissue.
4. To release attachment from tissue, turn the vacuum control stopcock OPEN to the atmosphere prior to removal of the stabilizer attachment.

WARNING

- Use of excessive negative vacuum pressure may cause tissue damage.
- Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.

STORAGE

The Titan Stabilizer Attachments storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

This page is intentionally left blank.

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Terumo® is a registered trademark of TERUMO CORPORATION
Titan™ is a trademark of Terumo Cardiovascular Systems Corporation
TERUMO AUSTRALIA PTY LTD
Terumo Corp.